

Model 20000 Storage Cabinet

Securely Store Your Endoscopes
for Optimal Readiness

- ▶ Air purge with programmable timer and HEPA filter allows for air movement through the endoscope's internal channels
- ▶ A locking mechanism provides peace of mind knowing your endoscopes are safe
- ▶ Holds up to eight endoscopes
- ▶ Endoscopes hang straight not coiled
- ▶ Endoscopes are held securely in place by a specially designed bracket
- ▶ Easy to clean and maintain
- ▶ Made of durable powder coated steel




The Model 20000 Storage Cabinet provides an environment for endoscope storage that is free of water-absorbent materials (e.g. laminated wood). A gasket surrounds the door opening. When the door is closed, the gasket is compressed providing a secured enclosure.



Model 20000

Storage Cabinet

Dimensions



WASSENBURG
medical

Wassenburg Medical Inc • 144 Railroad Drive, Ivyland, PA 18974
Info@wassenburgmedical.com • 215-364-1477
wassenburgmedicalinc.com

AD.0023, Revision 02212024